Intragastric Balloon
Lex-Bal



The Intragastric Balloon reduces food intake and thus leads to weight loss.

The foundation of this procedure is that the presence of a balloon inflated with liquid inside the stomach, which creates a feeling of fullness, reduces food intake and thus leads to weight loss.

The balloon, made of silicone rubber, a biocompatible-proven material, is inserted endoscopically through the mouth, while the patient is under sedation, and is filled with physiological saline up to a volume ranging from 400 to 600 cm3, according to the kind of patient being treated.

The Balloon must be removed after 6 months and may be replaced by a new one, if deemed necessary by the physician.

Advantages:

- This is an outpatient procedure and on average it takes 15-30 minutes



Models

- 400-600 cm3

LEXEL S.R.L. - Luis Saenz Peña 1937 I CP1135 CABA, Argentina
Phone. +54 011 4305-5617 / 5165 / 4304-5617 I Free line: 0800-444-0424

info@lexelmedical.com / www.lexelmedical.com/